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Three-in-one: Aldehyde-, alcohol- and phenol-free cleaning, disinfection and removal of grease from medical devices (acc. to MDD) and surfaces (acc. to BPD)

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- Gentle on materials
- Rapid and widespread effectiveness (incl. TbB)
- Listed by the German Society for Hygiene and Microbiology (DGHM)
- High yield
- Pleasant scent
- Removes grease and protein
- Easy to dose
- Non-corrosive
- Good wetting properties

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Disinfectant concentrate for the wiping of surfaces (walls, sanitary areas, fixtures and clothing); universally applicable, also for use in retirement homes and hospitals (acc. to BPD). Also suitable for the cleaning and disinfection of medical devices (acc. to MDD) such as handpieces and contra-angles, prostheses and silicone replicas. Ideally suited for surfaces that react sensitively to alcohol (e.g. acrylic glass). Effective against bacteria (incl. Tb), fungi and viruses (HBV/HIV). Listed by the German Society for Hygiene and Microbiology (DGHM).

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english

## Efficacy spectrum:

- bactericidal (incl. Tb)
- fungicidal
- effective against HBV/HIV

## Composition:

100 g Floracin® N CONCENTRATE contains:  
7.5g Bis(3-aminopropyl) dodecylamine  
surfactants  
excipients  
water

## Application:

Dilute the concentrate with the prescribed amount of tap water and apply the resulting solution to the surface to be disinfected by wiping. Ensure complete and thorough wetting!

Floracin® N CONCENTRATE is applied at the following concentrations:

### 3% concentration:

Add water to 30 ml concentrate to bring the volume to a total of 1 litre  
Inactivates bacteria and fungi in accordance with DGHM: exposure time 5 min.

### 2% concentration:

Add water to 20 ml concentrate to bring the volume to a total of 1 litre  
Inactivates bacteria and fungi in accordance with DGHM: exposure time 15 min.  
Inactivates HBV/HIV: exposure time 30 min. (in cases of higher organic contamination)

### 1.5% concentration:

Add water to 15 ml concentrate to bring the volume to a total of 1 litre  
Inactivates bacteria and fungi in accordance with DGHM: exposure time 30 min.

### 1% concentration:

Add water to 10 ml concentrate to bring the volume to a total of 1 litre  
Inactivates bacteria and fungi in accordance with DGHM: exposure time 60 min.  
Inactivates HBV/HIV: exposure time 30 min. (in cases of low organic contamination)

### 0.5% concentration:

Add water to 5 ml concentrate to bring the volume to a total of 1 litre  
Inactivates bacteria and fungi in accordance with DGHM: exposure time 240 min.  
Inactivates TbB (Mycobact. terrae): exposure time 60 min.

### 0.25% concentration:

Add water to 2.5 ml concentrate to bring the volume to a total of 1 litre  
Inactivates salmonellae:  
exposure time 5 min.

## Special advices:

Not suitable for use on parquet and laminate flooring or on fabric covers!

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FAVODENT Karl Huber GmbH, Greschbachstr. 17, 76229 Karlsruhe