

Guttapercha



Range of application:

Guttapercha is a temporary sealing compound for cavities with low bite pressure and that contain no liquid pharmaceuticals. When warmed, the guttapercha sticks can be easily molded, especially for pushing back the gums and for temporarily separating teeth.

Guttapercha has a high seal density and does not irritate dental pulp.

Temporary sealing compound in stick form

- Easy to mold
- High seal density
- No irritation of dental pulp
- Easy to remove

Composition:

100 g Guttapercha contain:

- 16.0 g guttapercha
- 81.0 g zinc oxide
- 2.4 g paraffi
- 1.2 g bees wax

Color: white

Unit	Art.-Nr.
32 g Package	501
192 g Package	502

On request we would be pleased to submit tailored conditions.

+49-721-62686-0
info@favodent.com